

HPLC Safety Set



GOES WITH
ALL **HPLC**
SYSTEMS!



One single Set for Solvent Supply ...



- With air valve
- Blocks solvent vapors
- Avoids underpressure
- Catches dust
- Keeps solvents clean
- Protects your health



... and Waste!

- With exhaust filter
- Blocks solvent vapors
- Avoids overpressure
- Protects your health
- Indicator shows filter status



Find more solutions at www.scatt-europe.com

Order Information

Find more solutions at www.scat-europe.com

CLEVER!

107 337 HPLC Safety Set



FOR USE WITH
ANY **HPLC**
SYSTEM!

Including:

107 910
SafetyCap III
GL 45

4x



117 010
Air valve

4x



501 113
Solvent
bottle
1 liter
GL 45
(glass)

4x



160 501
Blind plug
for capillary
connections

12x



117 808
Tube
connector
6.4-9 mm

1x



108 921
Safety
Waste Cap
GL 45

1x



610 535
Exhaust
filter
size M

1x



107 951
Container
5 liters
GL 45
PE-HD

1x



107 062
Fitting
Ø 2.3 mm
grey

10x



117 816
Tube
connector
6-8 mm

2x



107 063
Fitting
Ø 3.2 mm
blue

4x



107 064
Fitting
Ø 3.2 mm
red

4x



107 065
Fitting
Ø 3.2 mm
yellow

6x



107 061
Fitting
Ø 1.6 mm
green


10x



CERTAINLY THE
BEST SYSTEM
FOR HPLC LABS!

Discover the System!

Find more than 600 safety solutions in our catalogue.
Download for free at www.scat-europe.com
or ask your lab supplier for more information!





www.scat-europe.com

SCAT Europe GmbH • Opelstraße 3 • 64546 Mörfelden-Walldorf, Germany
Phone: +49-(0) 6105 - 305 586 - 0 • Fax: +49-(0) 6105 - 305 586 - 99

