



MS NOISE ENCLOSURES FOR BENCHTOP INSTRUMENTS (USB4 model)

Product Overview

This USB4 model reduces the typical noise of an operating ultrasonic unit to a tolerable level (noise reduction by approx. - 20dB(A)).

It allows operation of ultrasonic units in sensitive working areas, such as in laboratories, medical practices, manufacturing...

Vapours are carried off by an exhaust fan to protect the enclosure against an accumulation of heat / humidity.

Compatible with all major manufacturers such as Elma, Branson, ...

Key features :

- Transparent doors offer viewing convenience while limiting noise to an acceptable level.
- Easy double access to the instrument (from the front and also from the top. Top access is equipped with 2 gas springs for fast and frequent access)
- Decibel lowered by 20 dB(A)
- Opening for ultrasonic power supply and pipes placed at the back of the enclosure.
- Vapours are carried off by an exhaust fan to protect the enclosure against an accumulation of heat or humidity.
- Easy cleaning.





MS NOISE ENCLOSURES FOR BENCHTOP INSTRUMENTS (USB4 model)



Detailed Technical specifications :

- Environment friendly noise absorbing material, oil and water resistant.
- The noise reduction is in conformity with the European directive related to low voltage equipments (2006/95/EC) and electromagnetic compatibility (2004/108/EC).
- Compatible with 230V/50HZ or 110V/60HZ. delivered with a set of European, US, Japan, and UK plugs
- Power supply Csa, Ce and Ccc
- Electronic components are all compliant with RoHS directives (Restriction of Hazardous Substances Directive).

<i>External Dimensions</i>	Enclosure USB4 model (cm)	Enclosure USB4 model (inch)
Depth	Custom	Custom
Width	Custom	Custom
Height	Custom	Custom

Summary

Description	Part-number
Noise Reduction Enclosure (with 2x convenient access)	USB4CUSTxxxx

MS NOISE EUROPE & ASIA - MS NOISE - Building 4 ZAC des Vauguilletes - 1, bd des Noyers Pompons - F-89100 SENS – France (contact@msnoise.com)
MS NOISE NORTH AMERICA - MS NOISE - Suite 410 - 185, Alewife Brook Parkway - Cambridge, MA 02138 – USA (contact_US@msnoise.com)