

FORTUNA® POLYFIX® Bottle Top Dispenser



 Made
in
Germany

For all
chemicals
compatible
with PP!

101.072

101.071

101.070


Poulten & Graf
SUPERIOR LABORATORY PRODUCTS

FORTUNA®



Top Facts about FORTUNA® POLYFIX®

- dispense with German quality
- 3 POLYFIX® types for all chemicals compatible with PP
- visible dispensing procedure
- easy cleaning because of easy dismantling
- easy-to-use and solid volume setting system
- valve block made of polypropylene with FIX adapter
- comes with one PTFE intake tube and two PP adaptors (GL 40 + 45)



	Volume	grad.	Art.-No.
CLASSIC TYPE: FORTUNA® POLYFIX® with glass piston			
for all chemicals (incl. solvents), that do not attack polypropylene (PP) and do not crystallize	2-10 ml	0.2	101.070-38
	10-50 ml	1.0	101.070-48

ROBUST TYPE: FORTUNA® POLYFIX® with coated PTFE-piston			
for weak acids and alkaline solutions that do not attack polypropylene (PP)	2-10 ml	0.2	101.071-38
	10-50 ml	1.0	101.071-48

SHADY TYPE: FORTUNA® POLYFIX® with amber glass cylinder & glass piston			
for all light sensitive liquids that do not attack polypropylene (PP)	2-10 ml	0.2	101.072-38
	10-50 ml	1.0	101.072-48

Poulsen & Graf GmbH

Karl-Carstens-Straße 10
97877 Wertheim, Germany

www.poulsen-graf.de

sales@poulsen-graf.com

t + 49 9342 9229-0

f + 49 9342 9229-80

FORTUNA®