

FORTUNA UNIVERSAL[®] OPTIFIX[®] HF Bottle Top Dispenser



 Made
in
Germany

For almost all liquids -
including hydrofluoric
acid!

Poulten & Graf
SUPERIOR LABORATORY PRODUCTS

FORTUNA[®]



Top Facts about FORTUNA UNIVERSAL® OPTIFIX® HF:



- dispense liquids with German technology
 - valve system **without metal spring** - no contamination of dosing liquid
 - **suitable for analysis of micronutrients**
 - easy cleaning because of easy dismantling
 - easy-to-use and solid volume setting system
 - cylinder made of massive PTFE for extreme tightness
 - PTFE-encased piston with a sealing lip at its end for better tightness
-
- comes with one PTFE intake tube and two PTFE adaptors (GL 40 + 45)

Volume	grad.	Accuracy	Coefficient of Variation	Art.-No.
1 - 5 ml	0.2 ml	$\leq \pm 0,6 \%$	$\leq \pm 0.2 \%$	101.106-33
2 - 10 ml	0.2 ml	$\leq \pm 0,6 \%$	$\leq \pm 0.2 \%$	101.106-38
6 - 30 ml	1.0 ml	$\leq \pm 0,6 \%$	$\leq \pm 0.2 \%$	101.106-44
10 - 50 ml	1.0 ml	$\leq \pm 0,6 \%$	$\leq \pm 0.2 \%$	101.106-48