

Questron Sample Cart

Ease and Convenience for Storage and Transport of Samples



- Dedicated trolley for storing and transporting samples before and after digestion
- Can hold up to 840 samples placed in sample tray
- Provision to store up to 20 empty sample trays
- Sample distributed in 4 compartment accessible from two sides
- Ceramic tiled top provide convenient way to record samples in particular compartment
- Acid proof and corrosion free construction
- Heavy duty and high quality castor wheel for easy mobility.
- Inclined shelf to lock sample trays in place during movement
- Ideal for labs performing moderate to high volume of acid digestions
- Dimension Width x Depth x Height 71 x 71 x 127 cm (28 x 28 x 50 in)
- Custom size to accommodate your vials available

Questron provides total solution to your sample preparation needs
Microwave Digestion • Microwave Ashing • Hot Block Digestion
Automated Hot Block • Automated Workstation



Questron Technologies Corp.

6725 Millcreek Drive, Unit 7, Mississauga, Ontario L5N 5V3 CANADA
Tel: 905.363.1223 — Fax: 905.363.1227 — info@qtechcorp.com — www.qtechcorp.com